

Shell locking elbow joint passive device

6-01-EJ1




6-01-EJ1-CA

6-01-EJ1-DKBW

6-01-EJ1 Shell locking elbow joint passive device

- With passive elbow lock, and upper arm rotation joint



Product Code	Length	Proximal Circumference	Color	Wrist outer circumference	Wrist inner circumference
6-01-EJ1-CA	approx. 287mm	200mm	Caucasian	approx. 54mm	approx. 49mm
6-01-EJ1-DKBW	approx. 287mm	200mm	Dark brown	approx. 54mm	approx. 49mm

Cosmetic prostheses for forearm

6-01-EJ2



6-01-EJ2


6-01-EJ2 Cosmetic prostheses for forearm

- For above-elbow residual limb, with passive elbow lock, upper arm and forearm rotation, separate friction

Product Code	Side
6-01-EJ2-L	Left
6-01-EJ2-R	Right

Cosmetic prostheses for shoulder joint

6-01-SJ1



6-02-SJ1

6-01-SJ1 Cosmetic prostheses for shoulder joint

- For shoulder disarticulation, with passive elbow lock

Product Code	Side
6-01-SJ1-L	Left
6-01-SJ1-R	Right

Mechanical hand - voluntary opening (with PVCIH)

6-02-MH1



6-02-MH1

6-02-MH1 Mechanical hand voluntary opening (with PVCIH)

- With perlon cable through the back of the hand (outer pull)
- With threaded stud and inner hand

Product Code	Side	Size	Threaded stud
6-02-MH1-L-M	Left	M (7 1/4")	M12 x 1.5
6-02-MH1-R-M	Right	M (7 1/4")	M12 x 1.5
6-02-MH1-L-L	Left	L (7 3/4")	M12 x 1.5
6-02-MH1-R-L	Right	L (7 3/4")	M12 x 1.5
6-02-MH1-L-XL	Left	XL (8" - 8 1/4")	M12 x 1.5
6-02-MH1-R-XL	Right	XL (8" - 8 1/4")	M12 x 1.5

Pediatric

Feet & Adaptors

Knees & Hip Joints

Suspension Systems

Cosmetic, Tools & Materials

Orthotics & Rehabilitation