

## Eclipse Silicone Suspension Sleeve

3-01-ECSS



3-01-ECSS



### 3-01-ECSS Eclipse Silicone Suspension Sleeve

The Eclipse Silicone Suspension Sleeve provides extreme durability and comfort for your patient. Imagine not having to replace suspension sleeves every 3 months. A unique combination of materials and smart design and created a superior suspension sleeve that offers high performance and secure fit.

- Reinforced silicone and strategically located Parylene® coating provides maximum durability
- Dual sleeve design with non-slip inner strip of surface provides a maximum airtight seal which eliminates slippage
- This revolutionary patent pending design helps keep the sleeve in place even during the most strenuous activity
- Depending on activity, the outer Eclipse sleeve easily folds down over itself to reduce the length
- The Parylene® coating does not inhibit knee flexion which allows for excellent performance and comfort
- The Eclipse offers a platinum cured medical grade silicone with a Parylene® surface which reduces surface friction between the sleeve and skin helping in the prevention of skin shear for fragile skin and blister reduction
- The smooth and silky coating allows for clothing to naturally hang and not bunch up
- Warranty: 120 day warranty against defect and workmanship

**Measure circumference 15cm (6in) above mid-patella and choose the correct size**

Product Code	Color	Set contains	Size
<b>3-01-ECSS-S-BK</b>	Black	Long (outer)+Short (inner)	Small (41): 41cm – 49cm
<b>3-01-ECSS-M-BK</b>	Black	Long (outer)+Short (inner)	Medium (49): 49cm – 55cm
<b>3-01-ECSS-L-BK</b>	Black	Long (outer)+Short (inner)	Large (55): 55cm – 62cm
<b>3-01-ECSS-S-TN</b>	Tan	Long (outer)+Short (inner)	Small (41): 41cm – 49cm
<b>3-01-ECSS-M-TN</b>	Tan	Long (outer)+Short (inner)	Medium (49): 49cm – 55cm
<b>3-01-ECSS-L-TN</b>	Tan	Long (outer)+Short (inner)	Large (55): 55cm – 62cm